

# Schedule of Accreditation



Organisation Name Ireland East Hospital Group Our Lady's Hospital Navan  
 Trading As  
 INAB Reg No 215MT  
 Contact Name Ray O'Hare  
 Address Pathology Laboratory, Our Ladys Hospital, Navan, Meath  
 Contact Phone No 046 90 78571  
 Email ray.ohare@hse.ie  
 Website  
 Accreditation Standard EN ISO 15189  
 Standard Version 2012  
 Date of award of accreditation 17/02/2009  
 Scope Classification Blood Transfusion Science  
 Services available to the public<sup>1</sup>

<sup>1</sup> Refer to document on interpreting INAB Scopes of Accreditation

Sites from which accredited services are delivered		
(the detail of the accredited services delivered at each site are on the Scope of Accreditation)		
	Name	Address
1	Head Office	Pathology Laboratory, Navan, Meath

# Scope of Accreditation

## Head Office

### Blood Transfusion Science

Category: A

Medical pathology field - Test	Test/assay	Specimen Type	Equipment/Technique	Method (CE/Non-CE/In house developed/based on standard method)	Range of measurement	Std. ref & SOP
1020 Transfusion science - .01 Blood grouping including ABO, Rh(D) and other antigens by manual methods	Blood Grouping (ABO & RhD Typing)	Plasma	Tube method and automated column agglutination technology	Standard method - manual agglutination in tubes using CE marked reagents	Not applicable	LM-BT-0002 LM-BT-0003 LP-BT-0044
1020 Transfusion science - .03 Blood group antibody screen	Antibody Screening		Automated and manual method column agglutination technology	Standard method - manual Indirect Antiglobulin Test using CE marked Gel Cards	Not applicable	LM-BT-0004 LP-BT-0044
1020 Transfusion science - .04 Identification of blood group antibodies	Antibody Identification		Automated and manual method column agglutination technology	Standard method - manual Indirect Antiglobulin Test using CE marked Gel Cards & Enzyme Gel Cards.	Not applicable	LM-BT-0006 LP-BT-0044

1020 Transfusion science - .05 Cross match compatible donor units	Cross-Match Testing		Gel Technology	Standard method - manual Indirect Antiglobulin Test using CE marked Gel Cards	Not applicable	LM-BT-0015 Crossmatching using Gel Technology
1020 Transfusion science - .06 Red cell phenotyping	Antigen Typing	Red Blood Cells	Automated and manual column agglutination technology , tube method	Standard method - manual test using CE marked Gel Cards	Not applicable	LM-BT-0013 LM-BT-0014 LP-BT-0044
1020 Transfusion science - .09 Direct antiglobulin test	Direct Antiglobulin Test		Automated and manual method column agglutination technology	Standard method - manual Direct Antiglobulin Test using CE marked Gel Cards	Not applicable	LM-BT-0017 LP-BT-0044
1020 Transfusion science - .99 Miscellaneous tests	Electronic Issue	N/A	Electronic Issue on LIS	Electronic Issue on LIS	Not applicable	LP-BT-0045 Electronic Issue of Red Cells

*The hospital blood bank has been assessed and is competent to comply with Articles 14 and 15 of the EU Directive 2002/98/EC (S.I. 360/2005 and S.I. 547/2006)*