

# Schedule of Accreditation



Organisation Name	HSE National Drug Treatment Centre
Trading As	
INAB Reg No	169T
Contact Name	Maura Kehoe
Address	McCarthy Centre, 30-31 Pearse Street, Dublin
Contact Phone No	01 6488 645
Email	mkehoe@dtcb.ie
Website	
Accreditation Standard	ISO 17025 T
Date Initially Awarded	13/12/2005
Scope Classification	Chemical testing
Services available to the public <sup>1</sup>	Yes

<sup>1</sup> Refer to document on interpreting INAB Scopes of Accreditation

Sites from which accredited services are delivered	
(the detail of the accredited services delivered at each site are on the Scope of Accreditation)	
Name	Address
1 HSE National Drug Treatment Centre	McCarthy Centre, 30-31 Pearse Street, Dublin, D2

# Scope of Accreditation

## HSE National Drug Treatment Centre

### Chemical Testing

Category: A

Chemistry Field - Tests	Test name	Analyte	Range of measurement	Matrix	Equipment/technique	Standard reference/SOP
752 Chemical residue testing - .01 Drugs and drug metabolites	6-Monoacetyl Morphine	6-Monoacetyl Morphine	10 ng/ml	Urine	Beckman Coulter AU5800/Immunoassay	TP8
			8ng/ml	Oral Fluid	Beckman Coulter AU5800/Immunoassay	TP48
	Amphetamine	Amphetamine	50ng/ml	Oral Fluid	Beckman Coulter AU5800/Immunoassay	TP48
	Amphetamine/Ecstasy	Amphetamine/Ecstasy	1000ng/ml	Urine	Beckman Coulter AU5800/Immunoassay	TP8
	Benzodiazepines	Benzodiazepines	200ng/ml	Urine	Beckman Coulter AU5800/Immunoassay	TP8
			20ng/ml	Oral Fluid	Beckman Coulter AU5800/Immunoassay	TP48
	Buprenorphine	Buprenorphine	5ng/ml	Urine	Beckman Coulter AU5800/Immunoassay	TP8
	Cannabis	Cannabis	4ng/ml	Oral Fluid	Beckman Coulter AU5800/Immunoassay	TP48
	Cannabis Metabolites	Cannabis Metabolites	50ng/ml	Urine	Beckman Coulter AU5800/Immunoassay	TP8
	Cocaine	Cocaine	300ng/ml	Urine	Beckman Coulter AU5800/Immunoassay	TP8
	Cocaine		20ng/ml	Oral Fluid	Beckman Coulter AU5800/Immunoassay	TP48
	Creatinine	Creatinine	2ng/ml - 20ng/ml	Urine	Beckman Coulter AU5800/Immunoassay	TP8
	EDDP	EDDP	100ng/ml	Urine	Beckman Coulter AU5800/Immunoassay	TP8
	Ethyl Alcohol	Ethyl Alcohol	30mg/dl	Urine	Beckman Coulter AU5800/Immunoassay	TP8
	Ethyl Glucuronide	Ethyl Glucuronide	500ng/ml	Urine	Beckman Coulter AU5800/Immunoassay	TP8
	Methadone	Methadone	50ng/ml	Oral Fluid	Beckman Coulter AU5800/Immunoassay	TP48
	Methamphetamine	Methamphetamine	50ng/ml	Oral Fluid	Beckman Coulter AU5800/Immunoassay	TP48
	Opiates	Opiates	40ng/ml	Oral Fluid	Beckman Coulter AU5800/Immunoassay	TP48
			300ng/ml	Urine	Beckman Coulter AU5800/Immunoassay	TP8
	pH - Detect Test	pH - Detect Test	3 -11 pH units	Urine	Beckman Coulter AU5800/Immunoassay	TP8

	Serum Methadone	Methadone	0 - 500ng/ml	Serum	Randox Methadone ELISA Kit with BioTek LLx50 & BioTek ELx808	TP26
--	-----------------	-----------	--------------	-------	---	------