

Schedule of Accreditation



Organisation Name	Ireland East Hospital Group Our Lady's Hospital Navan
Trading As	
INAB Reg No	215MT
Contact Name	Ray O'Hare
Address	Pathology Laboratory, Our Ladys Hospital, Navan, Meath
Contact Phone No	046 90 78571
Email	ray.ohare@hse.ie
Website	
Accreditation Standard	EN ISO 15189
Standard Version	2012
Date of award of accreditation	17/02/2009
Scope Classification	Microbiology and virology
Scope Classification	Blood Transfusion Science
Services available to the public ¹	

¹ Refer to document on interpreting INAB Scopes of Accreditation

Sites from which accredited services are delivered		
(the detail of the accredited services delivered at each site are on the Scope of Accreditation)		
	Name	Address
1	Head Office	Pathology Laboratory, Navan, Meath

Scope of Accreditation

Head Office

Blood Transfusion Science

Category: A

Medical pathology field - Test	Test/assay	Specimen Type	Equipment/Technique	Method (CE/Non-CE/In house developed/based on standard method)	Range of measurement	Std. ref & SOP
1020 Transfusion science - .01 Blood grouping including ABO, Rh(D) and other antigens by manual methods	Blood Grouping (ABO & RhD Typing)	Plasma	Tube method and automated column agglutination technology	Standard method - manual agglutination in tubes using CE marked reagents	Not applicable	LM-BT-0002 LM-BT-0003 LP-BT-0044
1020 Transfusion science - .03 Blood group antibody screen	Antibody Screening		Automated and manual method column agglutination technology	Standard method - manual Indirect Antiglobulin Test using CE marked Gel Cards	Not applicable	LM-BT-0004 LP-BT-0044
1020 Transfusion science - .04 Identification of blood group antibodies	Antibody Identification		Automated and manual method column agglutination technology	Standard method - manual Indirect Antiglobulin Test using CE marked Gel Cards & Enzyme Gel Cards.	Not applicable	LM-BT-0006 LP-BT-0044
1020 Transfusion science - .05 Cross match	Cross-Match Testing		Gel Technology	Standard method - manual Indirect	Not applicable	LM-BT-0015 Crossmatching using

compatible donor units				Antiglobulin Test using CE marked Gel Cards		Gel Technology
1020 Transfusion science - .06 Red cell phenotyping	Antigen Typing	Red Blood Cells	Automated and manual column agglutination technology , tube method	Standard method - manual test using CE marked Gel Cards	Not applicable	LM-BT-0013 LM-BT-0014 LP-BT-0044
1020 Transfusion science - .09 Direct antiglobulin test	Direct Antiglobulin Test		Automated and manual method column agglutination technology	Standard method - manual Direct Antiglobulin Test using CE marked Gel Cards	Not applicable	LM-BT-0017 LP-BT-0044
1020 Transfusion science - .99 Miscellaneous tests	Electronic Issue	N/A	Electronic Issue on LIS	Electronic Issue on LIS	Not applicable	LP-BT-0045 Electronic Issue of Red Cells
<i>The hospital blood bank has been assessed and is competent to comply with Articles 14 and 15 of the EU Directive 2002/98/EC (S.I. 360/2005 and S.I. 547/2006)</i>						

Head Office

Microbiology and Virology

Category: A

Medical pathology field - Test	Test/assay	Specimen Type	Equipment/Technique	Method (CE/Non-CE/In house developed/based on standard method)	Range of measurement	Std. ref & SOP
1015 Detection and/or identification of bacterial, parasite, fungal and viral nucleic acids - .01 Nucleic acid probe hybridization, CE marked commercial systems	Automated in vitro diagnostic test for qualitative detection of Influenza A, Influenza B, and RSV	Nasopharangeal Swab	Cepheid GeneXpert	CE Marked	Detected/Not Detected	LP-MIC-0068 Detection of SAR-CoV-2, Flu A & B and RSV on GeneXpert
1015 Detection and/or identification of bacterial, parasite, fungal and viral nucleic acids - .03 Nucleic acid amplification tests, CE marked commercial systems	Sars-CoV-2	Nasophargeal Swab	GeneXpert	CE	Detected/Not Detected	LP-MIC-0068