

# Schedule of Accreditation



Organisation Name	The National Maternity Hospital
Trading As	
INAB Reg No	240MT
Contact Name	Damian Lally
Address	Department of Pathology and Laboratory Medicine, Holles Street, Dublin, D2
Contact Phone No	01 6373313
Email	dlally@nmh.ie
Website	
Accreditation Standard	EN ISO 15189
Standard Version	2022
Date of award of accreditation	15/09/2009
Scope Classification	Microbiology and Virology
Scope Classification	Blood Transfusion Science
Scope Classification	Haematology
Scope Classification	Histopathology and Cytopathology
Scope Classification	Chemical Pathology
Services available to the public <sup>1</sup>	Yes

<sup>1</sup> Refer to document on interpreting INAB Scopes of Accreditation

Sites from which accredited services are delivered		
(the detail of the accredited services delivered at each site are on the Scope of Accreditation)		
	Name	Address
1	Head Office	Department of Pathology and Laboratory Medicine, The National maternity Hospital, Holles Street, Dublin, Dublin, Ireland, 2

# Scope of Accreditation

## Head Office

### Blood Transfusion Science

Category: A

Medical pathology field - Test	Test/assay	Specimen Type	Equipment/Technique	Method (CE/Non-CE/In house developed/based on standard method)	Range of measurement	Std. ref & SOP
1020 Transfusion science - .01 Blood grouping including ABO, Rh(D) and other antigens by manual methods	ABO and Rh D Blood Grouping	Red Cells & Plasma	Manual	Manual	N/A	LP-BTR-ABORHT
1020 Transfusion science - .02 Blood grouping including ABO, Rh(D) and other antigens by automated methods			IH500	Automated	N/A	PP-CS-BT-30
1020 Transfusion science - .03 Blood group antibody screen	Antibody Screen	Plasma	IH500	Automated	N/A	PP-CS-BT-30
			Manual	Manual	N/A	LP-BTR-ABORHT

1020 Transfusion science - .04 Identification of blood group antibodies	Antibody Identification			Manual	N/A	LP-BTR-ABINVEST
1020 Transfusion science - .05 Cross match compatible donor units	Compatibility Testing	Plasma & Red Cells	1.Manual 2. Electronic Issue	1.Manual 2.Electronic Issue	N/A	LP-BTR-XMATCH, PP-CS-BT-28
1020 Transfusion science - .06 Red cell phenotyping	Antigen Typing	Red cells	IH500	Automated	N/A	PP-CS-BT-30
			Manual	Manual	N/A	LP-BTR-AT
1020 Transfusion science - .08 Antibody titration	Antibody Titre	Plasma and Red cells	IH500	Automated	N/A	PP-CS-BT-5
			Manual	Manual	N/A	PP-CS-BT-5
1020 Transfusion science - .09 Direct antiglobulin test	Direct Coombs Test	Red cells	IH500	Automated	N/A	PP-CS-BT-30
			Manual	Manual	N/A	LP-BTR-ABORHT

*The hospital blood bank has been assessed and is competent to comply with Articles 14 and 15 of the EU Directive 2002/98/EC (S.I. 360/2005 and S.I. 547/2006) OR (if applicable)*

## Head Office

### Chemical Pathology

Category: A

Medical pathology field - Test	Test/assay	Specimen Type	Technique	Equipment/Range of Measurement	Method (CE/Non-CE/In house developed/based on standard method)	Std. Ref & SOP
1061 Clinical Chemistry - .01 Analytes in general use in cardiac, liver function, lipid, renal and other profiles and metabolic studies	24 Hour Urinary Protein **1,2,3,4	Urine	Calculation	LIS	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
	Albumin **1,2,3,4	Plasma / Serum	Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
	Alkaline Phosphatase **1,2,3,4		Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
	ALT **1,2,3,4		Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
	Amylase **1,2,3,4		Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
	AST **1,2,3,4		Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure

					Analyser
Bilirubin-Direct **1,2,3,4		Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
Bilirubin-Total **1,2,3,4		Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
Calcium **1,2,3,4	Plasma	Arsenazo III	Beckman Coulter AU680	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
Chloride **1,2,3,4	Plasma / Serum / Urine	Indirect ISE	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
CK **1,2,3,4	Plasma / Serum	Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
Corrected Calcium**1,2,3,4		Calculation	LIS	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
Creatinine **1,2,3,4	Plasma / Serum / Urine	Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
Creatinine Clearance **1,2,3,4	Urine	Calculation	LIS	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
Glucose **1,2,3,4	CSF	Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the

					Roche Cobas Pure Analyser
	Fluoride	Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
IL-6 **1,2,3,4	Serum	Immunoassay	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
Inorganic Phosphate **1,2,3,4	Plasma	Phosphomolybdate Complex	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
LDH **1,2,3,4		LDH (L-P) IFCC	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
Magnesium **1,2,3,4	Plasma / Serum	Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
Osmolality **1,3,4	Plasma/Urine	Freezing-point depression	Advanced Micro Osmometer 3320	CE	PP-CS-BIO-34 Operation of the Advanced Micro-Osmometer Model 3320
Potassium **1,2,3,4	Plasma / Serum / Urine	Indirect ISE	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
Protein **1,2,3,4	CSF	Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser

	Urine	Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
Protein: Creatinine Ratio **1,2,3,4		Calculation	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser, PP-CS- BIO-42 Management of Reported Calculations in Biochemistry
Sodium **1,2,3,4	Plasma / Serum / Urine	Indirect ISE	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
Total Bile Acids **1,2,3,4	Plasma / Serum	Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
Total Protein **1,2,3,4		Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
Triglyceride **1,2,3,4		Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
Urate **1,2,3,4		Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
Urea **1,2,3,4	Plasma	GLDH, Kinetic Assay	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser

1061 Clinical Chemistry - .02 Proteins, quantitative analysis	CRP **1,2,3,4	Plasma / Serum	Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
1061 Clinical Chemistry - .10 Drugs for therapeutic monitoring	Gentamicin **1,2,3,4		Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
	Vancomycin **1,2,3,4		Electrochemiluminescence	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
1061 Clinical Chemistry - .20 Hormones	Anti-Mullerian Hormone **1,2,3,4	Serum	Immunoassay	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
	Free T4 **1,2,3,4	Serum / Plasma	Immunoassay	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
	HCG **1,2,3,4	Plasma	Immunoassay	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
	Oestradiol **1,2,3,4	Serum	Immunoassay	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
	PIGF **1,2,3,4	Serum	Immunoassay	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
	Procalcitonin **1,2,3,4		Immunoassay	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser

	Progesterone **1,2,3,4	Plasma	Immunoassay	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
	sFLT-1 **1,2,3,4	Serum	Immunoassay	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
	sFLT-1: PIGF ratio **1,2,3,4	Serum	Electrochemiluminescence immunoassay / Calculation	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser / PP-CS- BIO-42 Management of Reported Calculations in Biochemistry
	TSH **1,2,3,4	Serum / Plasma	Immunoassay	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
1061 Clinical Chemistry - .32 Maternal markers for foetal wellbeing	free $\beta$ -hCG **1,2,3,4	Serum	Immunoassay	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
	PAPP-A **1,2,3,4	Serum	Immunoassay	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
1061 Clinical Chemistry - .40 Iron studies	Ferritin **1,2,3,4	Serum	Immunoassay	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser
1061 Clinical Chemistry - .50 Protein and peptide tumour markers	CA-125 **1,2,3,4	Serum	Immunoassay	Roche Cobas Pure	CE	PP-CS-BIO-48 Operation of the Roche Cobas Pure Analyser

## Head Office

### Haematology

Category: A

Medical pathology field - Test	Test/Assay	Specimen Type	Technique	Range of Measurement/Equipment	Method (CE/Non-CE/In house developed/based on standard method)	Std. Ref & SOP
1030 Haematology - .01 Blood counts	Haematocrit **1,2,3,4	Whole Blood	Hydro Dynamic Focusing (DC Detection), flow cytometry	Sysmex XN-2000	CE	PP-CS-HAE-42 Operation of Sysmex FBC Analyser
	Haemoglobin **1,2,3,4		SLS-Haemoglobin Method	Sysmex XN-2000	CE	PP-CS-HAE-42 Operation of Sysmex FBC Analyser
	Mean Cell Haemoglobin **1,2,3,4		Hydro Dynamic Focusing (DC Detection), flow cytometry	Sysmex XN-2000	CE	PP-CS-HAE-42 Operation of Sysmex FBC Analyser
	Mean Cell Haemoglobin Concentration **1,2,3,4		Hydro Dynamic Focusing (DC Detection), flow cytometry	Sysmex XN-2000	CE	PP-CS-HAE-42 Operation of Sysmex FBC Analyser
	Mean Cell Volume **1,2,3,4		Hydro Dynamic Focusing (DC Detection), flow cytometry	Sysmex XN-2000	CE	PP-CS-HAE-42 Operation of Sysmex FBC Analyser
	Platelet Count **1,2,3,4		Hydro Dynamic Focusing (DC Detection), flow cytometry. Impedance / Optical	Sysmex XN-2000	CE	PP-CS-HAE-42 Operation of Sysmex FBC Analyser

	PLT-F **1,2,3,4		Hydro Dynamic Focusing (DC Detection), Flow Cytometry	Sysmex XN-2000	CE	PP-CS-HAE-42 Operation of Sysmex FBC Analyser
	Red Cell Count **1,2,3,4		Hydro Dynamic Focusing (DC Detection), flow cytometry	Sysmex XN-2000	CE	PP-CS-HAE-42 Operation of Sysmex FBC Analyser
	White Cell Count **1,2,3,4		Hydro Dynamic Focusing (DC Detection), flow cytometry	Sysmex XN-2000	CE	PP-CS-HAE-42 Operation of Sysmex FBC Analyser
1030 Haematology - .02 Visual examination of blood films	Basophils **1,2,3,4	Blood	Staining / Microscopy	Hematek stainer / Olympus Microscope	CE	PP-CS-HAE-17
	Eosinophils **1,2,3,4		Staining / Microscopy	Hematek stainer / Olympus Microscope	CE	PP-CS-HAE-17
	Lymphocytes **1,2,3,4		Staining / Microscopy	Hematek stainer / Olympus Microscope	CE	PP-CS-HAE-17
	Monocytes **1,2,3,4		Staining / Microscopy	Hematek stainer / Olympus Microscope	CE	PP-CS-HAE-17
	Neutrophils **1,2,3,4		Staining / Microscopy	Hematek stainer / Olympus Microscope	CE	PP-CS-HAE-17
	Platelet Morphology **1,2,3,4		Staining / Microscopy	Hematek stainer / Olympus Microscope	CE	PP-CS-HAE-17
	Red Cell Morphology **1,2,3,4		Staining / Microscopy	Hematek stainer / Olympus Microscope	CE	PP-CS-HAE-17
1030 Haematology - .05 Automated differential leucocyte counts	Basophil Count **1,2,3,4	Whole Blood	Hydro Dynamic Focusing (DC Detection), flow cytometry	Sysmex XN-2000	CE	PP-CS-HAE-42 Operation of Sysmex FBC Analyser
	Eosinophil Count **1,2,3,4		Hydro Dynamic Focusing (DC Detection), flow	Sysmex XN-2000	CE	PP-CS-HAE-42 Operation of

			cytometry			Sysmex FBC Analyser
	Lymphocyte Count **1,2,3,4		Hydro Dynamic Focusing (DC Detection), flow cytometry	Sysmex XN-2000	CE	PP-CS-HAE-42 Operation of Sysmex FBC Analyser
	Monocyte Count **1,2,3,4		Hydro Dynamic Focusing (DC Detection), flow cytometry	Sysmex XN-2000	CE	PP-CS-HAE-42 Operation of Sysmex FBC Analyser
	Neutrophil Count **1,2,3,4		Hydro Dynamic Focusing (DC Detection), flow cytometry	Sysmex XN-2000	CE	PP-CS-HAE-42 Operation of Sysmex FBC Analyser
	n-RBC's **1,3,4		Hydro Dynamic Focusing (DC Detection), flow cytometry	Sysmex XN-2000	CE	PP-CS-HAE-42 Operation of Sysmex FBC Analyser
1030 Haematology - .06 Automated reticulocyte counts	Reticulocyte Count **1,3,4		Hydro Dynamic Focusing (DC Detection), flow cytometry	Sysmex XN-2000	CE	PP-CS-HAE-42 Operation of Sysmex FBC Analyser
1030 Haematology - .30 Tests for haemoglobin variants and thalassaemia	Sickle Cell Screen **3,4	Blood	Qualitative Solubility Test	Manual	CE	PP-CS-HAE-13
1030 Haematology - .31 Tests for foetal Hb (Kleihauer)	Estimation of Foetal Maternal Haemorrhage - Kleihauer Test **1,2,4		Semi - Quantitative Acid Elution / Microscopy	Manual	CE	PP-CS-HAE-19
1030 Haematology - .41 General haemostasis related tests	Activated Partial Thromboplastin Time (APTT) **1,2,4		Coagulometric/Chromogenic	ACL Top 550	CE	PP-CS-HAE-43
	D-Dimer **1,2,4		Immuno-turbidimetric	ACL Top 550	CE	PP-CS-HAE-43

	Fibrinogen (Claus) **1,2,4	Coagulometric/Chromogenic	ACL Top 550	CE	PP-CS-HAE-43
	International Normalised Ratio (INR) **1,2,4	Coagulometric/Chromogenic	ACL Top 550	CE	PP-CS-HAE-43
	Prothrombin Time (PT) **1,2,4	Coagulometric/Chromogenic	ACL Top 550	CE	PP-CS-HAE-43
1030 Haematology - .57 Screening test for infectious mononucleosis	Detection of Infectious mononucleosis heterophile antibodies **3,4	Chromatographic immunoassay for the qualitative detection of IM heterophile antibodies	Manual	CE	PP-CS-HAE-8

## Head Office

### Histopathology and Cytopathology

Category: A

Medical pathology field - Test	Test/assay	Specimen Type	Equipment/Technique	Method (CE/Non-CE/In house developed/based on standard method)	Range of measurement	Std. ref & SOP
1051 Histopathology - .01 Processing fixed specimens for Histopathological testing	Cut Up, Tissue Processing, Embedding, Microtomy, Automated and Manual Haemotoxylin & Eosin staining, Coverslipping, Semi-Automated and Automated Immunohistochemical Tissue Staining	Human Tissue	Histocore Pearl / VIP5JR / Roche HE600 Autostainer / Ventana / Roche Benchmark Ultra	CE	N/A	LP-AP-TISSUEPROC LP-AP-TISSUEPROC PP-CS-AP-44 LP-AP-EMBEDDING LP-AP-Microtomy HM325 LP-AP-AUTOMICROTOMY LP-AP-AUTOSTAIN PP-CS-AP-45 PP-CS-AP-64 PP-CS-AP-65
1051 Histopathology - .03 Histochemistry	Special stains - Gram Stain Oil Red O PAS PAS-D Perls Prussian Blue Manual H+E		Manual	CE	N/A	LP-AP-GRAM, LP-AP-OILREDO, LP-AP-PAS, LPAP-PERLS, LP-AP-Manual H+E
1051 Histopathology - .05 Histological interpretation- paediatric pathology	Diagnostic Interpretation and reporting of Human tissue Specialist areas excluded:				N/A	PP-CS-AP-13, PP-CS-AP-10, PP-CS-AP-33, PP-CS-AP-11, PP-CS-AP-34, PP-CS-AP-37, PP-CS-AP-12, PP-CS-AP-24, PP-CS-AP-35,

	Medical Renal Pathology Paediatric Pathology Neuropathology Breast screening Pathology specimens Specialist areas included: Perinatal Pathology (confined to Placental examination & Products of conception)
1051 Histopathology - .09 Immunohistochemistry	Actin - SMA
	AE1/AE3
	BCL 2
	CA 125
	Calretinin
	CAM5.2
	CD 10
	CD 20
	CD 3
	CD 34
CD 56	

			PP-CS-AP-32, PP-CS-AP-9, PP-CS-AP-14, PP-CS-AP-2, PP-CS-AP-28, PP-CS-AP-47
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche	CE	N/A	PP-CS-AP-45, PP-CS-

CD138
CD31
CD68
CEA
Chromogranin A
CK 5/6
CK20
CK7
CMV
D240
Desmin
EMA
GFAP
HCGp
HMB45
HPL

Benchmark Ultra			AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	Non-CE/In House Developed	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Ventana, Roche	CE	N/A	PP-CS-AP-45, PP-CS-

	Benchmark Ultra			AP-65
Inhibin	Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
K167/P16 Dual Stain	Ventana Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45
Ki-67	Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Oestrogen Receptor	Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
P 53	Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
P16	Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
P57	Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Parvovirus	Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
PLAP	Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Progesterone Receptor	Ventana Benchmark Ultra	CE	N/A	PP-CS-AP-45
S100	Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Synaptophysin	Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
Vimentin	Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65
WT1	Ventana, Roche Benchmark Ultra	CE	N/A	PP-CS-AP-45, PP-CS-AP-65

## Head Office

### Microbiology and Virology

Category: A

Medical pathology field - Test	Test/assay	Specimen Type	Equipment/Technique	Method (CE/Non-CE/In house developed/based on standard method)	Range of measurement	Std. ref & SOP
1011 Macroscopic examination and description	Blood Culture ** 2,3,4	Blood	BacT Alert	CE	N/A	PP-CS-MIC-2
			BacT Alert	CE	N/A	PP-CS-MIC-2
	CSF Culture ** 2,3,4	CSF	Culture and sensitivity	CE	N/A	PP-CS-MIC-12
	Sputum Culture ** 2,3,4	Endotracheal aspirates, sputum	Culture and sensitivity	CE	N/A	PP-CS-MIC-37
	Tip Cultures ** 2,3,4	UVC Tip, UAC Tip, Long Line Tip, PICC Line Tip, CVP Tip, Urinary Catheter Tip, Central Line Tip, Endotracheal Tube Tip, PEG Line Tip, PAC Tip, Arterial Line Tip, Chest Drain Tip	Culture and sensitivity	CE	N/A	PP-CS-MIC-19, PP-CS-MIC-25
1012 Preparation of films on glass slides followed by microscopic examination with or without fixation and staining with dyes as required - .01 Microscopic	Blood Culture ** 4	Blood	BacT Alert, Previ Colour, microscopy	CE	N/A	PP-CS-MIC-2, PP-CS-MIC-24, PP-CS-MIC-9

examination for general bacteriology purposes (including enumeration and description of human cells)						
	CSF Culture ** 4	CSF	Microscopy, Previ Colour	CE	N/A	PP-CS-MIC-12, PP-CS-MIC-24
	Ear, Nose, Throat and Oral Culture ** 4	nasal swab, nasal secretions, sinus aspirate, antral washout, mouth swab, tongue swab, ear swab, mastoid swab, throat swab	Microscopy, Previ Colour	CE	N/A	PP-CS-MIC-15, PP-CS-MIC-24
	Eye Culture ** 4	Eye swab, canalicular pus, corneal scraping, corneal swab, contact lens, Vitreous\Aqueous tap, orbital cellulitis swab	Microscopy, Previ Colour	CE	N/A	PP-CS-MIC-20, PP-CS-MIC-24
	Films for Microscopic Examination ** 4	Films for microscopic examination	Microscopy, Previ Colour	CE	N/A	PP-CS-MIC-24, PP-CS-MIC-9
	Fluid Cultures ** 4	Amniotic fluid, Ascitic Fluid, Bursa Fluid, Pericardial Fluid, Peritoneal fluid, Pleural fluid, Synovial fluid	Previ Colour, microscopy	CE	N/A	PP-CS-MIC-38, PP-CS-MIC-24
	Genital Culture ** 4	High vaginal swab (HVS), vaginal discharge, vulval swab, labial swab, endocervical swab, penile swab, low vaginal swab, urethral	Microscopy, Previ Colour	CE	N/A	PP-CS-MIC-23, PP-CS-MIC-24, PP-CS-MIC-9

		swab, cervical swab, semen, buttock swab, groin swab.				
	Sputum Culture ** 4	Endotracheal aspirates, sputum	Microscopy, Previ Colour	CE	N/A	PP-CS-MIC-19, PP-CS-MIC-24, PP-CS-MIC-37
	Urine Culture ** 4	Bag Urine, MSU, CSU, Supra-Pubic, Clean Catch	Microscopy	CE	N/A	PP-CS-MIC-42
	Wound Culture, Abscesses and Deep Seated Infections, Skin Infections ** 4	Exudates, Pus, Wound Swabs, Abscesses	Microscopy, Previ Colour	CE	N/A	PP-CS-MIC-44, PP-CS-MIC-24
1012 Preparation of films on glass slides followed by microscopic examination with or without fixation and staining with dyes as required - .03 Microscopic examination for fungi	Blood Culture ** 4	Blood	BacT Alert, microscopy	CE	N/A	PP-CS-MIC-2, PP-CS-MIC-34
	CSF Culture ** 4	CSF	Microscopy	CE	N/A	PP-CS-MIC-12, PP-CS-MIC-34
	Ear, Nose, Throat and Oral Culture ** 4	Nasal swab, nasal secretions, sinus aspirate, antral washout, mouth swab, tongue swab, ear swab, mastoid swab, throat swab	Microscopy	CE	N/A	PP-CS-MIC-15, PP-CS-MIC-34
	Eye Culture ** 4	Eye swab, canalicular pus, corneal scraping, corneal swab, , contact lens, Vitreous/Aqueous tap, orbital cellulitis swab	Microscopy	CE	N/A	PP-CS-MIC-20, PP-CS-MIC-34

	Films for microscopic examination ** 4	Films for microscopic examination	Microscopy	CE	N/A	PP-CS-MIC-34, PP-CS-MIC-57
	Fluid Cultures ** 4	Amniotic fluid, Ascitic Fluid, Bursa Fluid, Pericardial Fluid, Peritoneal fluid, Pleural fluid, Synovial fluid	Microscopy	CE	N/A	PP-CS-MIC-38, PP-CS-MIC-34
	Genital Culture ** 4	High vaginal swab (HVS), vaginal discharge, vulval swab, labial swab, endocervical swab, penile swab, low vaginal swab, urethral swab, cervical swab, semen, buttock swab, groin swab.	Microscopy	CE	N/A	PP-CS-MIC-23, PP-CS-MIC-34
	Limited Identification of Fungal Isolates ** 4	Body fluids, excretions, materials, swabs and environmental and infection control swabs	Microscopy	CE	N/A	PP-CS-MIC-18, PP-CS-MIC-34
	Sputum Culture ** 4	Endotracheal aspirates, sputum	Microscopy	CE	N/A	PP-CS-MIC-19, PP-CS-MIC-34, PP-CS-MIC-37
	Wound Culture, Abscesses and Deep Seated Infections, Skin Infections ** 4	Exudates, Pus, Wound Swabs, Abscesses, Skin Swab	Microscopy	CE	N/A	PP-CS-MIC-44, PP-CS-MIC-34
1013 Culture of organisms in liquid or agar based culture media with visual or instrument monitoring for growth -	Blood Culture ** 2,3,4	Blood	BacT Alert, culture and sensitivity	CE	N/A	PP-CS-MIC-2, PP-CS-MIC-28

.01 Culture of general bacteria						
	CSF Culture ** 2,3,4	CSF	Culture and sensitivity	CE	N/A	PP-CS-MIC-12, PP-CS-MIC-28
	Ear, Nose, Throat and Oral Culture ** 2,3,4	Ear swab, mastoid swab, nasal swab, nasal secretions, sinus aspirate, antral washout, throat swab, mouth swab, tongue swab	Culture and sensitivity	CE	N/A	PP-CS-MIC-15, PP-CS-MIC-28
	Expressed Breast Milk, Breast Swabs and Nipple Swabs ** 2,3,4	Expressed breast milk, breast swabs and nipple swabs	Culture and sensitivity	CE	N/A	PP-CS-MIC-17, PP-CS-MIC-28
	Eye Culture ** 2,3,4	Eye swab, canalicular pus, corneal scraping, corneal swab, orbital cellulitis swab, contact lens, Vitreous/Aqueous tap	Culture and sensitivity, BacT Alert	CE	N/A	PP-CS-MIC-20, PP-CS-MIC-2, PP-CS-MIC-28
	Fluid Cultures ** 2,3,4	Amniotic fluid, Acsitic Fluid, Bursa Fluid, Pericardial Fluid, Peritoneal fluid, Pleural fluid, Synovial fluid	Culture and sensitivity	CE	N/A	PP-CS-MIC-28, PP-CS-MIC-38
	Genital Culture ** 2,3,4	High vaginal swab (HVS), vaginal discharge, vulval swab, labial swab, endocervical swab, penile swab, low vaginal swab, urethral swab, cervical swab, semen, buttock swab,	Culture and sensitivity	CE	N/A	PP-CS-MIC-23, PP-CS-MIC-28

	groin swab.				
MRSA Screen ** 2,3,4	Nasal, throat, axilla, groin, perineum, umbilical, eyes and other sites depending on clinical information.	Culture and sensitivity	CE	N/A	PP-CS-MIC-29, PP-CS-MIC-28
Placental Swabs ** 2,3,4	Placental swabs	Culture and sensitivity	CE	N/A	PP-CS-MIC-33, PP-CS-MIC-28
PM Screening Swabs and IUD Swabs ** 2,3,4	PM screening swabs, PM CSF, IUD swabs	Culture and sensitivity	CE	N/A	PP-CS-MIC-6, PP-CS-MIC-28
Screening of Rectal Swabs for Carriage of VRE, CPE and Gentamicin Resistant Gram- Negative Bacilli ** 2,3,4	Faeces, rectal swab	Culture and sensitivity	CE	N/A	PP-CS-MIC-8, PP-CS-MIC-28
Sputum Culture ** 2,3,4	Endotracheal aspirates, sputum	Culture and sensitivity	CE	N/A	PP-CS-MIC-19, PP-CS-MIC-37, PP-CS-MIC-28
Tip Cultures ** 2,3,4	UVC Tip, UAC Tip, Long Line Tip, PICC Line Tip, CVP Tip, Urinary Catheter Tip, Central Line Tip, Endotracheal Tube Tip, PEG Line Tip, PAC Tip, Arterial Line Tip, Chest Drain Tip	Culture and sensitivity	CE	N/A	PP-CS-MIC-19, PP-CS-MIC-25, PP-CS-MIC-28
Urine Culture ** 2,3,4	Bag Urine, MSU, CSU, Supra-Pubic, Clean Catch	Microscopy, culture and sensitivity	CE	N/A	PP-CS-MIC-42, PP-CS-MIC-28
Wound Culture, Abscesses and Deep	Exudates, Pus, Wound Swabs,	Culture and sensitivity	CE	N/A	PP-CS-MIC-44, PP-CS-MIC-28

	Seated Infections, Skin Infections ** 2,3,4	Abscesses, Skin Swab				
1013 Culture of organisms in liquid or agar based culture media with visual or instrument monitoring for growth - .02 Culture of fungi	Blood Culture ** 2,3,4	Blood	BacT Alert, culture and sensitivity	CE	N/A	PP-CS-MIC-2, PP-CS-MIC-28
	CSF Culture ** 2,3,4	CSF	Culture and sensitivity	CE	N/A	PP-CS-MIC-12, PP-CS-MIC-28
	Ear, Nose, Throat and Oral Culture ** 4	Ear swab, mastoid swab, nasal swab, nasal secretions, sinus aspirate, antral washout, throat swab, mouth swab, tongue swab	Culture and sensitivity	CE	N/A	PP-CS-MIC-15, PP-CS-MIC-28
	Expressed Breast Milk, Breast Swabs and Nipple Swabs ** 4	Expressed breast milk, breast swabs and nipple swabs	Culture and sensitivity	CE	N/A	PP-CS-MIC-17, PP-CS-MIC-28
	Eye Culture ** 4	Eye swab, canalicular pus, corneal scraping, corneal swab, orbital cellulitis swab, contact lens, Vitreous/Aqueous tap	Culture and sensitivity, BacT Alert	CE	N/A	PP-CS-MIC-20, PP-CS-MIC-2, PP-CS-MIC-28
	Fluid Cultures ** 4	Amniotic fluid, Acsitic Fluid, Bursa Fluid, Pericardial Fluid, Peritoneal fluid, Pleural fluid, Synovial fluid	Culture and sensitivity	CE	N/A	PP-CS-MIC-28, PP-CS-MIC-38
	Genital Culture ** 4	High vaginal swab (HVS), vaginal	Culture and sensitivity	CE	N/A	PP-CS-MIC-23, PP-CS-MIC-28,

		discharge, vulval swab, labial swab, endocervical swab, penile swab, low vaginal swab, urethral swab, cervical swab, semen, buttock swab, groin swab.				
	Limited Identification of Fungal Isolates ** 4	Body fluids, excretions, materials, swabs and environmental and infection control swabs	Culture and sensitivity	CE	N/A	PP-CS-MIC-18, PP-CS-MIC-28
	Placental Swabs ** 4	Placental swabs	Culture and sensitivity	CE	N/A	PP-CS-MIC-33, PP-CS-MIC-28
	PM Screening Swabs and IUD Swabs ** 4	PM screening swabs, PM CSF, IUD swabs	Culture and sensitivity	CE	N/A	PP-CS-MIC-6, PP-CS-MIC-28
	Sputum Culture ** 4	Endotracheal aspirates, sputum	Culture and sensitivity	CE	N/A	PP-CS-MIC-19, PP-CS-MIC-37, PP-CS-MIC-28
	Tip Cultures ** 4	UVC Tip, UAC Tip, Long Line Tip, PICC Line Tip, CVP Tip, Urinary Catheter Tip, Central Line Tip, Endotracheal Tube Tip, PEG Line Tip, PAC Tip, Arterial Line Tip, Chest Drain Tip	Culture and sensitivity	CE	N/A	PP-CS-MIC-19, PP-CS-MIC-25, PP-CS-MIC-28
	Wound Culture, Abscesses and Deep Seated Infections, Skin Infections ** 4	Exudates, Pus, Wound Swabs, Abscesses, Skin Swab	Culture and sensitivity	CE	N/A	PP-CS-MIC-44, PP-CS-MIC-28
1015 Detection and/or identification of bacterial,	BioFire Blood Culture Identification 2	Blood	FilmArray Multiplex PCR System	CE	Detected or Not Detected	PP-CS-MIC-77

parasite, fungal and viral nucleic acids - .03 Nucleic acid amplification tests, CE marked commercial systems	(BCID2) Panel ** 2,3,4					
	BioFire Respiratory Panel ** 2,3,4	Nasal Swabs, NPA	FilmArray Multiplex PCR System	CE	Detected or Not Detected	PP-CS-MIC-77
	Detection of CPE using Carba R on the GeneXpert ** 2,3,4	Bacterial isolates	GeneXpert	CE	Detected or Not Detected	PP-CS-MIC-76
	Detection of Faecal pathogens on BioFire GI Panel **1,2,3,4	Faeces	FilmArray Multiplex PCR System	CE	Detected / Not Detected	PP-CS-MIC-77, PP-CS-MIC-21
	Detection of GBS on the GeneXpert ** 2,3,4	HVS/Rectal Swab	GeneXpert	CE	Detected or Not Detected	PP-CS-MIC-76
	Detection of MRSA-SA on the GeneXpert ** 2,3,4	Blood Cultures, bacterial isolates	GeneXpert	CE	Detected or Not Detected	PP-CS-MIC-76
	Detection of Sars-COV-2 on the Aus Dx ** 2,3,4	Nasal Swabs, NPA	Aus Dx	CE	Detected or Not Detected	PP-CS-MIC-80
	Detection of Sars-CoV-2 Influenza A, B and RSV PCR using Xpert Xpress CoV-2 / Flu / RSV on GeneXpert **1,2,3,4	Nasal Swabs, NPA	GeneXpert	CE	Detected or Not Detected	PP-CS-MIC-76
	Xpert MRSA NxG Assay ** 2,3,4	Bacterial Isolates	GeneXpert	CE	Detected or Not Detected	PP-CS-MIC-76
1016 Identification of cultured bacteria and fungi using non-nucleic acid based techniques -	Blood Culture ** 2,3,4	Blood	BacT Alert, Previ Colour, culture and sensitivity, VITEK 2 and manual methods	CE	N/A	PP-CS-MIC-2, PP-CS-MIC-24, PP-CS-MIC-57, PP-CS-MIC-75

.01 Biochemical methods, CE marked commercial systems						
	CSF Culture ** 4	CSF	Culture and sensitivity, Previ Colour, VITEK 2 and manual	CE	N/A	PP-CS-MIC-12, PP-CS-MIC-24, PP-CS-MIC-57, PP-CS-MIC-75
	Ear, Nose, Throat and Oral Culture ** 2,3,4	Ear swab, mastoid swab, nasal swab, nasal secretions, sinus aspirate, antral washout, throat swab, mouth swab, tongue swab	Culture and sensitivity, Previ Colour, VITEK 2 and manual	CE	N/A	PP-CS-MIC-15, PP-CS-MIC-24, PP-CS-MIC-57, PP-CS-MIC-75
	Expressed Breast Milk, Breast Swabs and Nipple Swabs ** 2,3,4	Expressed breast milk, breast swabs and nipple swabs	Previ Colour, Culture and sensitivity, VITEK 2 and manual	CE	N/A	PP-CS-MIC-17, PP-CS-MIC-24, PP-CS-MIC-57, PP-CS-MIC-75
	Eye Culture ** 2,3,4	Eye swab, canaliculitis, corneal scraping, corneal swab, orbital cellulitis swab, contact lens, Vitreous/Aqueous tap	Culture and sensitivity, Previ Colour, VITEK 2 and manual	CE	N/A	PP-CS-MIC-20, PP-CS-MIC-24, PP-CS-MIC-57, PP-CS-MIC-75
	Fluid Cultures ** 2,3,4	Amniotic fluid, Acsitic Fluid, Bursa Fluid, Pericardial Fluid, Peritoneal fluid, Pleural fluid, Synovial fluid	Previ Colour, culture and sensitivity, VITEK 2 and manual	CE	N/A	PP-CS-MIC-38, PP-CS-MIC-24, PP-CS-MIC-57, PP-CS-MIC-75
	Genital Culture ** 2,3,4	High vaginal swab (HVS), vaginal discharge, vulval swab, labial swab, endocervical swab, penile swab, low	Culture and sensitivity, Previ Colour, VITEK 2 and manual	CE	N/A	PP-CS-MIC-23, PP-CS-MIC-24, PP-CS-MIC-57, PP-CS-MIC-75

	vaginal swab, urethral swab, cervical swab, semen, buttock swab, groin swab.				
MRSA Screen ** 2,3,4	Nasal, throat, axilla, groin, perineum, umbilical, eyes and other sites depending on clinical information.	Culture and sensitivity, Previ Colour, VITEK 2 and manual	CE	N/A	PP-CS-MIC-29, PP-CS-MIC-24, PP-CS-MIC-57, PP-CS-MIC-75
Placental Swabs ** 2,3,4	Placental swabs	Previ Colour, culture and sensitivity, VITEK 2 and manual	CE	N/A	PP-CS-MIC-33, PP-CS-MIC-24, PP-CS-MIC-57, PP-CS-MIC-75
PM Screening Swabs and IUD Swabs ** 2,3,4	PM screening swabs, PM CSF, IUD swabs	Previ Colour, culture and sensitivity, VITEK 2 and manual	CE	N/A	PP-CS-MIC-6, PP-CS-MIC-24, PP-CS-MIC-57, PP-CS-MIC-75
Screening of Rectal Swabs for Carriage of VRE, CPE and Gentamicin Resistant Gram- Negative Bacilli ** 2,3,4	Faeces, rectal swab	Culture and sensitivity, Previ Colour, VITEK 2 and manual	CE	N/A	PP-CS-MIC-8, PP-CS-MIC-24, PP-CS-MIC-57, PP-CS-MIC-75, PP-CS-MIC-48
Sputum Culture ** 2,3,4	Endotracheal aspirates, sputum	Previ colour, culture and sensitivity, VITEK 2 and manual	CE	N/A	PP-CS-MIC-19, PP-CS-MIC-37, PP-CS-MIC-24, PP-CS-MIC-57, PP-CS-MIC-75,
Tip Cultures ** 2,3,4	UVC Tip, UAC Tip, Long Line Tip, PICC Line Tip, CVP Tip, Urinary Catheter Tip, Central Line Tip, Endotracheal Tube Tip, PEG Line Tip,	Culture and sensitivity, Previ Colour, VITEK 2 and manual	CE	N/A	PP-CS-MIC-19, PP-CS-MIC-25, PP-CS-MIC-24, PP-CS-MIC-57, PP-CS-MIC-75

		PAC Tip, Arterial Line Tip, Chest Drain Tip				
	Urine Culture ** 2,3,4	Bag Urine, MSU, CSU, Supra-Pubic, Clean Catch	Culture and sensitivity, Previ Colour, VITEK 2 and manual	CE	N/A	PP-CS-MIC-42, PP-CS-MIC-24, PP-CS-MIC-57, PP-CS-MIC-75
	Wound Culture, Abscesses and Deep Seated Infections, Skin Infections ** 2,3,4	Exudates, Pus, Wound Swabs, Abscesses, Skin Swab	Culture and sensitivity, Previ Colour, VITEK 2 and manual	CE	N/A	PP-CS-MIC-44, PP-CS-MIC-24, PP-CS-MIC-57, PP-CS-MIC-75
1016 Identification of cultured bacteria and fungi using non-nucleic acid based techniques - .03 Identification of fungi by microscopic morphology	Blood Culture ** 4	Blood	Microscopy, culture and sensitivity, manual	CE	N/A	PP-CS-MIC-2, PP-CS-MIC-34, PP-CS-MIC-57
	CSF Culture ** 2,3,4	CSF	Microscopy, culture and sensitivity and manual	CE	N/A	PP-CS-MIC-12, PP-CS-MIC-34, PP-CS-MIC-57
	Ear, Nose, Throat and Oral Culture ** 4	Ear swab, mastoid swab, nasal swab, nasal secretions, sinus aspirate, antral washout, throat swab, mouth swab, tongue swab	Culture and sensitivity, microscopy and manual	CE	N/A	PP-CS-MIC-15, PP-CS-MIC-34, PP-CS-MIC-57
	Expressed Breast Milk, Breast Swabs and Nipple Swabs ** 4	Expressed breast milk, breast swabs and nipple swabs	Microscopy, culture and sensitivity and manual	CE	N/A	PP-CS-MIC-17, PP-CS-MIC-34, PP-CS-MIC-57
	Eye Culture ** 4	Eye swab, canalicular pus, corneal scraping, corneal swab, orbital cellulitis swab,	Microscopy, culture and sensitivity and manual	CE	N/A	PP-CS-MIC-20, PP-CS-MIC-34, PP-CS-MIC-57

	contact lens, Vitreous\Aqueous tap				
Fluid Cultures ** 4	Amniotic fluid, Acsitic Fluid, Bursa Fluid, Pericardial Fluid, Peritoneal fluid, Pleural fluid, Synovial fluid	Microscopy, culture and sensitivity, manual	CE	N/A	PP-CS-MIC-38, PP-CS-MIC-34, PP-CS-MIC-57
Genital Culture ** 4	High vaginal swab (HVS), vaginal discharge, vulval swab, labial swab, endocervical swab, penile swab, low vaginal swab, urethral swab, cervical swab, semen, buttock swab, groin swab.	Microscopy, culture and sensitivity and manual	CE	N/A	PP-CS-MIC-23, PP-CS-MIC-34, PP-CS-MIC-57
Limited Identification of Fungal Isolates ** 4	Body fluids, excretions, materials, swabs and environmental and infection control swabs	Microscopy, culture and sensitivity and manual	CE	N/A	PP-CS-MIC-18, PP-CS-MIC-34, PP-CS-MIC-57
Placental Swabs ** 4	Placental swabs	Microscopy, culture and sensitivity and manual	CE	N/A	PP-CS-MIC-33, PP-CS-MIC-34, PP-CS-MIC-57
PM Screening Swabs and IUD Swabs ** 4	PM screening swabs, PM CSF, IUD swabs	Microscopy, culture and sensitivity and manual	CE	N/A	PP-CS-MIC-6, PP-CS-MIC-34, PP-CS-MIC-57
Sputum Culture ** 4	Endotracheal aspirates, sputum	Microscopy, culture and sensitivity and manual	CE	N/A	PP-CS-MIC-19, PP-CS-MIC-37, PP-CS-MIC-34,
Tip Cultures ** 4	UVC Tip, UAC Tip, Long Line Tip, PICC Line Tip, CVP	Microscopy, culture and sensitivity and manual	CE	N/A	PP-CS-MIC-19, PP-CS-MIC-25, PP-CS-MIC-34, PP-CS-

		Tip, Urinary Catheter Tip, Central Line Tip, Endotracheal Tube Tip, PEG Line Tip, PAC Tip, Arterial Line Tip, Chest Drain Tip				MIC-57
	Wound Culture, Abscesses and Deep Seated Infections, Skin Infections ** 4	Exudates, Pus, Wound Swabs, Abscesses, Skin Swab	Microscopy, culture and sensitivity and manual	CE	N/A	PP-CS-MIC-44, PP-CS-MIC-34, PP-CS-MIC-57
1016 Identification of cultured bacteria and fungi using non-nucleic acid based techniques - .04 Identification using MALDI-TOF Spectroscopy	Blood Culture **1,2,3,4	Blood	Vitek MALDITOF MS	CE	N/A	PP-CS-MIC-79
	CSF Culture **1,2,3,4	CSF	Vitek MALDITOF MS	CE	N/A	PP-CS-MIC-79
	Ear, Nose, Throat and Oral Culture **1,2,3,4	Ear swab, mastoid swab, nasal swab, nasal secretions, sinus aspirate, antral washout, throat swab, mouth swab, tongue swab	Vitek MALDITOF MS	CE	N/A	PP-CS-MIC-79
	Environmental and Infection Control Swabs **1,2,3,4	Environmental fluid specimens, swabs, settle plates, contact plates, air sampling plates	Vitek MALDITOF MS	CE	N/A	PP-CS-MIC-79
	Expressed Breast Milk, Breast Swabs and Nipple Swabs **1,2,3,4	Expressed breast milk, breast swabs and nipple swabs	Vitek MALDITOF MS	CE	N/A	PP-CS-MIC-79
	Eye Culture **1,2,3,4	Eye swab, canalicular pus,	Vitek MALDITOF MS	CE	N/A	PP-CS-MIC-79

	corneal scraping, corneal swab, orbital cellulitis swab, contact lens, Vitreous Aqueous tap				
Fluid Cultures **1,2,3,4	Amniotic fluid, Acsitic Fluid, Bursa Fluid, Pericardial Fluid, Peritoneal fluid, Pleural fluid, Synovial fluid	Vitek MALDITOF MS	CE	N/A	PP-CS-MIC-79
Genital Culture **1,2,3,4	High vaginal swab (HVS), vaginal discharge, vulval swab, labial swab, endocervical swab, penile swab, low vaginal swab, urethral swab, cervical swab, semen, buttock swab, groin swab.	Vitek MALDITOF MS	CE	N/A	PP-CS-MIC-79
MRSA Screen **1,2,3,4	Nasal, throat, axilla, groin, perineum, umbilical, eyes and other sites depending on clinical information.	Vitek MALDITOF MS	CE	N/A	PP-CS-MIC-79
Placental Swabs **1,2,3,4	Placental swabs	Vitek MALDITOF MS	CE	N/A	PP-CS-MIC-79
PM Screening Swabs and IUD Swabs **1,2,3,4	PM screening swabs, PM CSF, IUD swabs	Vitek MALDITOF MS	CE	N/A	PP-CS-MIC-79
Screening of Rectal	Faeces, rectal swab	Vitek MALDITOF MS	CE	N/A	PP-CS-MIC-79

	Swabs for Carriage of VRE, CPE, ESBL and Gentamicin Resistant Gram Negative Bacilli **1,2,3,4					
	Sputum Culture **1,2,3,4	Endotracheal aspirates, sputum	Vitek MALDITOF MS	CE	N/A	PP-CS-MIC-79
	Tip Cultures **1,2,3,4	UVC Tip, UAC Tip, Long Line Tip, PICC Line Tip, CVP Tip, Urinary Catheter Tip, Central Line Tip, Endotracheal Tube Tip, PEG Line Tip, PAC Tip, Arterial Line Tip, Chest Drain Tip	Vitek MALDITOF MS	CE	N/A	PP-CS-MIC-79
	Urine Culture **1,2,3,4	Bag Urine, MSU, CSU, Supra-Pubic, Clean Catch	Vitek MALDITOF MS	CE	N/A	PP-CS-MIC-79
	Wound Culture, Abscesses and Deep Seated Infections, Skin Infections **1,2,3,4	Exudates, Pus, Wound Swabs, Abscesses, Skin Swab	Vitek MALDITOF MS	CE	N/A	PP-CS-MIC-79
1017 Measurement of antimicrobial activity and application of clinical interpretive criteria to general bacteria (rapidly growing aerobes) - .05 Other categories of organism (as specified)	Blood Culture **1,2,3,4	Blood	Culture and sensitivity, VITEK 2 and manual.	CE	N/A	PP-CS-MIC-2, PP-CS-MIC-48, PP-CS-MIC-57, PP-CS-MIC-75
	CSF Culture **1,2,3,4	CSF	Culture and sensitivity,	CE	N/A	PP-CS-MIC-12, PP-

		VITEK 2 and manual			CS-MIC-48, PP-CS-MIC-57, PP-CS-MIC-75
Ear, Nose, Throat and Oral Culture **1,2,3,4	Ear swab, mastoid swab, nasal swab, nasal secretions, sinus aspirate, antral washout, throat swab, mouth swab, tongue swab	Culture and sensitivity	CE	N/A	PP-CS-MIC-15, PP-CS-MIC-48, PP-CS-MIC-57, PP-CS-MIC-75
Expressed Breast Milk, Breast Swabs and Nipple Swabs **1,2,3,4	Expressed breast milk, breast swabs and nipple swabs	Culture and sensitivity, VITEK 2 and manual	CE	N/A	PP-CS-MIC-17, PP-CS-MIC-48, PP-CS-MIC-57, PP-CS-MIC-75
Eye Culture **1,2,3,4	Eye swab, canalicular pus, corneal scraping, corneal swab, orbital cellulitis swab, contact lens, Vitreous/Aqueous tap	Culture and sensitivity	CE	N/A	PP-CS-MIC-20, PP-CS-MIC-48, PP-CS-MIC-57, PP-CS-MIC-75
Fluid Cultures **1,2,3,4	Amniotic fluid, Acsitic Fluid, Bursa Fluid, Pericardial Fluid, Peritoneal fluid, Pleural fluid, Synovial fluid	Culture and sensitivity, VITEK 2 and manual	CE	N/A	PP-CS-MIC-38, PP-CS-MIC-48, PP-CS-MIC-57, PP-CS-MIC-75
Genital Culture **1,2,3,4	High vaginal swab (HVS), vaginal discharge, vulval swab, labial swab, endocervical swab, penile swab, low vaginal swab, urethral swab, cervical swab, semen, buttock swab, groin swab.	Culture and sensitivity, VITEK 2 and manual	CE	N/A	PP-CS-MIC-23, PP-CS-MIC-48, PP-CS-MIC-57, PP-CS-MIC-75

MRSA Screen **1,2,3,4	Nasal, throat, axilla, groin, perineum, umbilical, eyes and other sites depending on clinical information.	Culture and sensitivity, VITEK 2 and manual	CE	N/A	PP-CS-MIC-29, PP-CS-MIC-48, PP-CS-MIC-57, PP-CS-MIC-75
Placental Swabs **1,2,3,4	Placental swabs	Culture and sensitivity, VITEK 2 and manual	CE	N/A	PP-CS-MIC-33, PP-CS-MIC-48, PP-CS-MIC-57, PP-CS-MIC-75
PM Screening Swabs and IUD Swabs **1,2,3,4	PM screening swabs, PM CSF, IUD swabs	Culture and sensitivity, VITEK 2 and manual	CE	N/A	PP-CS-MIC-6, PP-CS-MIC-48, PP-CS-MIC-57, PP-CS-MIC-75
Screening of Rectal Swabs for Carriage of VRE, CPE and Gentamicin Resistant Gram- Negative Bacilli **1,2,3,4	Faeces, rectal swab	Culture and sensitivity, VITEK 2 and manual	CE	N/A	PP-CS-MIC-8, PP-CS-MIC-48, PP-CS-MIC-57, PP-CS-MIC-75
Sputum Culture **1,2,3,4	Endotracheal aspirates, sputum	Culture and sensitivity	CE	N/A	PP-CS-MIC-19, PP-CS-MIC-37, PP-CS-MIC-48, PP-CS-MIC-57, PP-CS-MIC-75
Tip Cultures **1,2,3,4	UVC Tip, UAC Tip, Long Line Tip, PICC Line Tip, CVP Tip, Urinary Catheter Tip, Central Line Tip, Endotracheal Tube Tip, PEG Line Tip, PAC Tip, Arterial Line Tip, Chest Drain Tip	Culture and sensitivity, VITEK 2 and manual	CE	N/A	PP-CS-MIC-19, PP-CS-MIC-25, PP-CS-MIC-48, PP-CS-MIC-57, PP-CS-MIC-75
Urine Culture **1,2,3,4	Bag Urine, MSU, CSU, Supra-Pubic,	Culture and Sensitivity, Manual and VITEK 2	CE	N/A	PP-CS-MIC-42, PP-CS-MIC-48, PP-CS-

		Clean Catch				MIC-57, PP-CS-MIC-75
	Wound Culture, Abscesses and Deep Seated Infections, Skin Infections **1,2,3,4	Exudates, Pus, Wound Swabs, Abscesses, Skin Swab	Culture and sensitivity, VITEK 2 and manual	CE	N/A	PP-CS-MIC-44, PP-CS-MIC-48, PP-CS-MIC-57, PP-CS-MIC-75
1024 Preservation of microbial cultures	Blood Culture ** 4	Blood	Culture and sensitivity	CE	N/A	PP-CS-MIC-2, PP-CS-MIC-57
	CSF Culture ** 4	CSF	Culture and sensitivity	CE	N/A	PP-CS-MIC-12, PP-CS-MIC-57
	Ear, Nose, Throat and Oral Culture ** 4	Ear swab, mastoid swab, nasal swab, nasal secretions, sinus aspirate, antral washout, throat swab, mouth swab, tongue swab	Culture and sensitivity	CE	N/A	PP-CS-MIC-15, PP-CS-MIC-57
	Expressed Breast Milk, Breast Swabs and Nipple Swabs ** 4	Expressed breast milk, breast swabs and nipple swabs	Culture and sensitivity	CE	N/A	PP-CS-MIC-17, PP-CS-MIC-57
	Eye Culture ** 4	Eye swab, canalicular pus, corneal scraping, corneal swab, orbital cellulitis swab, contact lens, Vitreous/Aqueous tap	Culture and sensitivity	CE	N/A	PP-CS-MIC-20, PP-CS-MIC-57
	Fluid Cultures ** 4	Amniotic fluid, Acsitic Fluid, Bursa Fluid, Pericardial Fluid, Peritoneal fluid, Pleural fluid, Synovial fluid	Culture and sensitivity	CE	N/A	PP-CS-MIC-38, PP-CS-MIC-57
	Genital Culture ** 4	High vaginal swab (HVS), vaginal	Culture and sensitivity	CE	N/A	PP-CS-MIC-23, PP-CS-MIC-57

	discharge, vulval swab, labial swab, endocervical swab, penile swab, low vaginal swab, urethral swab, cervical swab, semen, buttock swab, groin swab.				
Limited Identification of Fungal Isolates ** 4	Body fluids, excretions, materials, swabs and environmental and infection control swabs	Culture and sensitivity	CE	N/A	PP-CS-MIC-18, PP-CS-MIC-34, PP-CS-MIC-57.
MRSA Screen ** 4	Nasal, throat, axilla, groin, perineum, umbilical, eyes and other sites depending on clinical information.	Culture and sensitivity	CE	N/A	PP-CS-MIC-29, PP-CS-MIC-57
Placental Swabs ** 4	Placental swabs	Culture and sensitivity	CE	N/A	PP-CS-MIC-33, PP-CS-MIC-57
PM Screening Swabs and IUD Swabs ** 4	PM screening swabs, PM CSF, IUD swabs	Culture and sensitivity	CE	N/A	PP-CS-MIC-6, PP-CS-MIC-57
Screening of Rectal Swabs for Carriage of VRE, CPE and Gentamicin Resistant Gram- Negative Bacilli ** 4	Faeces, rectal swab	Culture and sensitivity	CE	N/A	PP-CS-MIC-8, PP-CS-MIC-57
Sputum Culture ** 4	Endotracheal aspirates, sputum	Culture and sensitivity	CE	N/A	PP-CS-MIC-19, PP-CS-MIC-37
Tip Cultures ** 4	UVC Tip, UAC Tip, Long Line Tip, PICC Line Tip, CVP Tip ,Urinary Catheter	Culture and sensitivity	CE	N/A	PP-CS-MIC-19, PP-CS-MIC-25, PP-CS-MIC-57

		Tip, Central Line Tip, Endotracheal Tube Tip, PEG Line Tip, PAC Tip, Arterial Line Tip, Chest Drain Tip				
	Urine Culture ** 4	Bag Urine, MSU, CSU, Supra-Pubic, Clean Catch	Microscopy, culture and sensitivity	CE	N/A	PP-CS-MIC-42, PP-CS-MIC-57
	Wound Culture, Abscesses and Deep Seated Infections, Skin Infections ** 4	Exudates, Pus, Wound Swabs, Abscesses, Skin Swab	Culture and sensitivity	CE	N/A	PP-CS-MIC-44, PP-CS-MIC-57
1029 Miscellaneous - .99 Miscellaneous tests	Pregnancy test **1,2,3,4	Urine	Manual, Rapid response	CE	N/A	PP-CS-MIC-50

*The laboratory has been awarded flexible scope in the scope classifications as noted in the scope document and in accordance with the laboratory's approved and documented procedures.*

*Note 1 - Range may be extended for the test*

*Note 2 – New parameters/tests may be added*

*Note 3 – New matrices may be added*

*Note 4 – Changes to equipment/kits where the underlying methodology does not change*

*For further details please refer to the laboratory's 'List of flexible scope changes', available directly from the laboratory.*