

Schedule of Accreditation



Organisation Name	Rotunda Hospital
Trading As	
INAB Reg No	208MT
Contact Name	Susan Luke
Address	Department of Laboratory Medicine, Parnell Square, Dublin, D1
Contact Phone No	01 8171739
Email	sluke@rotunda.ie
Website	
Accreditation Standard	EN ISO 15189
Standard Version	2012
Date of award of accreditation	09/12/2008
Scope Classification	Microbiology and virology
Scope Classification	Blood Transfusion Science
Scope Classification	Haematology
Scope Classification	Histopathology and cytopathology
Scope Classification	Chemical pathology

Services available to the public¹

¹ Refer to document on interpreting INAB Scopes of Accreditation

Sites from which accredited services are delivered		
(the detail of the accredited services delivered at each site are on the Scope of Accreditation)		
	Name	Address
1	Department of Laboratory Medicine	Department of Laboratory Medicine-Quality Office, Parnell Square, Dublin, D1

Scope of Accreditation

Department of Laboratory Medicine

Blood Transfusion Science

Category: A

Medical pathology field - Test	Test/assay	Specimen Type	Equipment/Technique	Method (CE/Non-CE/In house developed/based on standard method)	Range of measurement	Std. ref & SOP
1020 Transfusion science - .01 Blood grouping including ABO, Rh(D) and other antigens by manual methods	ABO and Rh D grouping of blood	Red Cells Plasma / Serum	Manual method.	Column Agglutination Technology (CAT) Microtitre Plate, Tube Technology	N/A	LM-BT-0002
1020 Transfusion science - .02 Blood grouping including ABO, Rh(D) and other antigens by automated methods	ABO and RhD grouping of blood	Red cells Plasma Serum	Ortho Vision	Column Agglutination Technology (CAT) CE Marked	N/A	LP-BT-0063
1020 Transfusion science - .03 Blood group antibody screen	Antibody Screening	Plasma / Serum	Automated Ortho Vision.	Column Agglutination Technology (CAT) , Tube Technology	N/A	LP-BT-0063
		Plasma and serum	Manual method	Column Agglutination Technology (CAT) CE Marked	N/A	LM-BT-0003 & LM-BT-0022

1020 Transfusion science - .04 Identification of blood group antibodies	Antibody Identification	Plasma / Serum	Manual/BioVueand Vision	Column Agglutination Technology (CAT)		LM-BT-0004
1020 Transfusion science - .05 Cross match compatible donor units	Compatibility Testing		Manual / Biorad BioVue	Column Agglutination Technology (CAT)		LM-BT-0022 LM-BT-0010
1020 Transfusion science - .06 Red cell phenotyping	Red Cell Phenotyping	Red Cells	Manual BioVue and Tube	Column Agglutination Technology (CAT) , Tube Technology		LM-BT-0006
1020 Transfusion science - .09 Direct antiglobulin test	Direct Antibody Test		Manual BioVue (CAT) and Tube	Column Agglutination Technology (CAT) , Tube Technology CE marked	N/A	LM-BT-0005
				Ortho Vision	Column Agglutination Technology (CAT) CE marked	N/A
1020 Transfusion science - .99 Miscellaneous tests	Antibody Titre	Plasma / Serum	Manual bioVue and Vision	Column Agglutination Technology (CAT)		LM-BT-0013
	FMH Estimation	EDTA red Cells	DX Flex flow cytometer	CE marked BRAD/AEVZ		LP-BT-0066, LM-BT-0019
The hospital blood bank has been assessed and is competent to comply with Articles 14 and 15 of the EU Directive 2002/98/EC (S.I. 360/2005 and S.I. 547/2006)						

Category: A

Medical pathology field - Test	Test/assay	Specimen Type	Technique	Equipment/Range of Measurement	Method (CE/Non-CE/In house developed/based on standard method)	Std. Ref & SOP
1061 Clinical Chemistry - .01 Analytes in general use in cardiac, liver function, lipid, renal and other profiles and metabolic studies	AFP ***1234	Serum	Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	Albumin **1234	Lithium Heparin	Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	ALP **1234		Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	ALT		Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	ALT **1234		Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Amylase **1234		Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	AST		Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	AST **1234		Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Bile Acids **1234		Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Calcium **1234		Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Chloride **24	Lithium Heparin / Urine	Indirect ISE	Roche Cobas Pure	CE Method	LP-BIO-0034
	CK **1234	Lithium Heparin	Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Corrected Calcium Calculation **1234		Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Creatinine **1234	Lithium Heparin / Urine	Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Creatinine Clearance Calculation	24hr Urine Collection	Calculated Test	Roche Cobas Pure	CE Method	LP-BIO-0034

	CRP **1234	Lithium Heparin	Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Direct Bilirubin **1234		Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Folate **1234	Serum	Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	Fructosamine/Corrected Fructosamine**1234	Lithium Heparin / serum	Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	GGT **1234	Lithium Heparin	Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Glucose **1234	Sodium Flouride/CSF	Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	IL6 **1234	Lithium Heparin	Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	Lactate **1234	Sodium Fluoride/CSF	Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	LDH **1234	Lithium Heparin	Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Magnesium **1234		Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Phosphate **1234		Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Potassium **1234	Lithium Heparin / Urine	Indirect ISE	Roche Cobas Pure	CE Method	LP-BIO-0034
	Protein Creatinine Ratio **1,2,3,4	Urine	Calculation	Roche Cobas Pure / LIMS	CE METHOD	LP-BIO-0020
	PTH **1234	Serum	Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	Sodium **1234	Lithium Heparin / Urine	Indirect ISE	Roche Cobas Pure	CE Method	LP-BIO-0034
	Total Bilirubin **1234	Lithium Heparin	Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Total Protein **1234	Lithium Heparin / Urine / CSF	Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Triglycerides **1234	Lithium Heparin	Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Urate **1234		Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Urea **1234		Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Urinary Protein	Urine	Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
1061 Clinical Chemistry - .10 Drugs for therapeutic monitoring	Gentamicin **1234	Lithium Heparin	Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034
	Vancomycin**1234		Electrochemiluminescence	Roche Cobas Pure	CE Method	LP-BIO-0034

1061 Clinical Chemistry - .20 Hormones	AMH **1234	Serum	Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	Bioactive Prolactin	Serum	Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	FAI**1,2,3,4	Serum	Calculated Test	Roche Cobas Pure	CE Method	LP-BIO-0034
	FSH **1234		Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	FT3		Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	FT4 **1234	Serum/Lithium Heparin	Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	HCG-B **1234		Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	LH **1234	Serum	Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	Oestradiol **1234		Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	Procalcitonin	Serum & Lithium Heparin	Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	Progesterone **1234	Serum	Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	Prolactin **1234		Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	SHBG **1234		Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	Testosterone **1234		Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	TSH **1234	Serum / Lithium Heparin	Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
1061 Clinical Chemistry - .40 Iron studies	Ferritin **1234	Serum / Lithium Heparin	Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
1061 Clinical Chemistry - .45 Vitamin B12 and folate	B12 **1234	Serum	Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	Folate **1234		Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
1061 Clinical Chemistry - .47 Vitamin assays	Vitamin D**1234		Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
1061 Clinical Chemistry - .50 Protein and peptide tumour markers	CA-125 **1234		Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	CA19.9		Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	CEA		Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034

	HE4		Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034
	ROMA Index Calculation		Calculated Test	Roche Cobas Pure	CE METHOD	LP-BIO-0034
1061 Clinical Chemistry - .61 Hb A1c	Haemoglobin A1C **1,2,3,4	Whole Blood EDTA	HPLC	HA-8180	CE METHOD	LP-BIO-0019
1061 Clinical Chemistry - .85 Anti-thyroglobulin antibodies	Anti-TPO **1234	Serum	Immunoassay	Roche Cobas Pure	CE Method	LP-BIO-0034

Category: B

Medical pathology field - Test	Test/assay	Specimen Type	Technique	Equipment/Range of Measurement	Method (CE/Non-CE/In house developed/based on standard method)	Std. Ref & SOP
1061 Clinical Chemistry - .01 Analytes in general use in cardiac, liver function, lipid, renal and other profiles and metabolic studies	Glucose **1234	Whole Blood	Enzymatic	Roche AccuChek Inform II	CE Method	LP-POC-0027
1061 Clinical Chemistry - .05 CO-oximetry	Bilirubin **2,4	Whole Blood - heparinised	Spectrophotometry	Radiometer ABL90 Flex 0-690umol/L	CE	LP-POC-0035
	CA15.3	Serum	Electrochemiluminescence Immunoassay	Roche Cobas Pure (1.5-300U/mL)	CE Method	LP-BIO-0034
	Hb **2,4	Whole Blood - heparinised	Spectrophotometry	Radiometer ABL90 Flex 0-270g/L	CE	LP-POC-0035
	MethHb**2,4,		Spectrophotometry	Radiometer ABL90 Flex 0-100%	CE	LP-POC-0035
1061 Clinical Chemistry - .06 Blood pH and gas tensions	Anion gap **2,4	Whole Blood-heparinised	Calculated	Radiometer ABL90 Flex	CE	LP-POC-0035
	Base Excess **2,4	Whole Blood - heparinised	Calculated	Radiometer ABL90 Flex	CE	LP-POC-0035
	HCO3- **2,4	Whole Blood - heparinised	Calculated	Radiometer ABL90 Flex	CE	LP-POC-0035
	pCO2 **2,4		Potentiometry	Radiometer ABL90 Flex 1.6-14.7kPa	CE	LP-POC-0035
	pH **2,4		Potentiometry	Radiometer ABL90 Flex 6.750-7.850	CE	LP-POC-0035
	pO2 **2,4		Optical pO2	Radiometer ABL90 Flex 1.33-73.3kPa	CE	LP-POC-0035

1061 Clinical Chemistry - .07 Other analytes performed on a blood gas analyser	Ca++ **2,4		Potentiometry (Direct ISE)	Radiometer ABL90 Flex 0.10- 2.70mmol/L	CE	LP-POC-0035
	Cl- **2,4		Potentiometry (Direct ISE)	Radiometer ABL90 Flex 70-160mmol/L	CE	LP-POC-0035
	Glucose **2,4		Amperometry	Radiometer ABL90 Flex 0-47 mmol/L	CE	LP-POC-0035
	K+ **2,4		Potentiometry (Direct ISE)	Radiometer ABL90 Flex 1.5-10.5mol/L	CE	LP-POC-0035
	Lactate **2,4		Amperometry	Radiometer ABL90 Flex 0-31	CE	LP-POC-0035
	Na+ **2,4		Potentiometry (Direct ISE)	Radiometer ABL90 Flex 94-190mmol/L	CE	LP-POC-0035
1061 Clinical Chemistry - .35 Assays of neonatal specimens	Glucose **1234	Whole Blood	Enzymatic	Roche AccuChek Inform II	CE Method	LP-POC-0027
1061 Clinical Chemistry - .76 Simple side tests for biochemical and immunological analytes	Glucose**1234		Roche AccuChek Inform II Poc Testing	Roche AccuChek Inform II	CE Enzymatic	LP-POC-0027 Operation of the Roche AccuChek Inform II Glucose Meter

Category: A

Medical pathology field - Test	Test/Assay	Specimen Type	Technique	Range of Measurement/Equipment	Method (CE/Non-CE/In house developed/based on standard method)	Std. Ref & SOP
1030 Haematology - .01 Blood counts	Full Blood Count including reticulocytes, white cell differential (automated) includes neutrophils, lymphocytes, eosinophils, monocytes and basophils. *1,2,3,4	Whole Blood, EDTA, Na Citrate, Thromboexact blood	Impedance counting, fluorescence analysis by flow cytometry and SLS haemoglobin measurement on automated cell counter.	Sysmex XN 2000	CE based on Standard method	LP-HAEM-0032
1030 Haematology - .02 Visual examination of blood films	Blood Film and white cell differential (manual) includes neutrophils, lymphocytes, monocytes, eosinophils, basophils and red cell morphology. **2,4	EDTA Blood	Microscopic examination with or without fixation and staining with dyes for enumeration and description of human cells	Aerospray stainer / microscopy	CE Marked	LM-HAEM-0007, 0008, LP-HAEM-0003, 0004, WI-HAEM-0003, LP-HAEM-0011 & 0006
1030 Haematology - .09 Examination of malarial parasites	Malaria Parasites Examination**1,2,4	Blood Film	Microscopic examination with or without fixation and staining with dyes for enumeration and description of parasites	Manual Microscopy	CE Marked	LM-HAEM-0009, 0010 & LP-HAEM-0006
	Rapid Diagnosis of Malaria *1,2,4	EDTA Blood	Immuno - chromatographic using monoclonal antibodies	Manual technique	Immunochromogenic using monoclonal antibodies	LM-HAEM-0010

1030 Haematology - .30 Tests for haemoglobin variants and thalassaemia	Adult and Paediatric Haemoglobinopathy screening**1,2,4	EDTA Blood / Clotted Blood	Capillary Electrophoresis	Sebia	CE Marked	LM-HAEM-0033 and -0034, LP-HAEM- 0028
1030 Haematology - .31 Tests for foetal Hb (Kleihauer)	Kleihauer (tests for foetal Hb)**1,2,4	EDTA Blood	Microscopic examination with or without fixation and staining with dyes for enumeration and description of human cells	Manual Microscopy	CE Marked	LM-HAEM-0021 & LP-HAEM-0004

Haematology

Category: B

Medical pathology field - Test	Test/Assay	Specimen Type	Technique	Range of Measurement/Equipment	Method (CE/Non-CE/In house developed/based on standard method)	Std. Ref & SOP
1030 Haematology - .01 Blood counts	Haemaglobin estimation**1,2,4	Whole Blood	Photometric	Hemocue 210DM	Photometric CE Marked	LP-POC-0009
1030 Haematology - .40 Limited haemostasis related tests	Lupus Anticoagulant	Sodium Citrate	Coagulometric/Chromogenic	dRVVTScreen: 16.0-240.0 seconds dRRVT confirm: 6.0-121.0 seconds	CE MARKED	LP-HAEM-0023, LM-HAEM-0031
1030 Haematology - .41 General haemostasis related tests	APTT/APTT50 **1,2,3,4		Coagulometric/Chromogenic	16.0-400 seconds / ACL TOP 550	CE MARKED	LP-HAEM-0023
	Clauss Fibrinogen **1,2,3,4		Coagulometric/Chromogenic	Fib-C 0.8-5.0 g/l, Fib-C L: 0.35-1.20 g/L, Fib-C H: 4.0-10.0 g/L / ACL TOP 550	CE MARKED	LP-HAEM-0023
	PT/PT50 **1,2,3,4		Coagulometric/Chromogenic	8.0-320.0 seconds/ACL TOP 550	CE MARKED	LP-HAEM-0023

Category: A

Medical pathology field - Test	Test/assay	Specimen Type	Equipment/Technique	Method (CE/Non-CE/In house developed/based on standard method)	Range of measurement	Std. ref & SOP
1051 Histopathology - .01 Processing fixed specimens for Histopathological testing	Automated Haematoxylin & Eosin **2,4	Human tissue/ Cell Block	DAKEWE DP360 Stainer & DAKEWE CS500 Coverslipper	CE Marked instrument	N/A	LM-HIST-0008 & LP-HIST-0079
	Cut-Up *2*4		Thermoscientific Printmate, digital recording system, Kern balance	CE Marked instrument	N/A	LP-HIST-0009, 0056,0069,0072
	Embedding *2*4		Leica Arcadia, embedding centre, heated forceps	CE MARKED	N/A	LP-HIST-0030 & LP-HIST-0072
	Microtomy *2*4		microtomes, waterbaths, incubators	CE Marked instrument	N/A	LP-HIST-0071 & 0072
	Tissue Processing **2,4		Sakura Tissue Tek VIP6 (E2) & Sakura Tissue Tek VIP 6 (AI)	CE Marked instrument	N/A	LP-HIST-0046 & 0081
1051 Histopathology - .03 Histochemistry	Special Stains: Manual H&E	Human tissue / Cell Block	Manual, DAKEWE DP360 Stainer and DP360 CS500 Coverslipper	Based on Standard Method	N/A	LM-HIST-0008 & LP-HIST-0079
	Special Stains: Papanicolaou		Manual, DAKEWE DP360 Stainer and DP360 CS500 Coverslipper	Based on Standard Method	N/A	LM-HIST-0008 & LP-HIST-0079
1051 Histopathology - .06 Histological interpretation-perinatal pathology	Diagnostic interpretation and reporting of human tissue specimens -		Microscope	CE Marked	N/A	LP-HIST-0032, LP-HIST-0073

	specialist areas included: perinatal pathology (confined to histological examination of placental tissue and gravid uteri). Specialist areas excluded: medical renal pathology, paediatric pathology, neuropathology & breast screening pathology)					
1051 Histopathology - .09 Immunohistochemistry	CD10 **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE Marked Instrument	N/A	LP-HIST-0080, LP-HIST-0079
	CD138 **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE Marked	N/A	LP-HIST-0080, LP-HIST-0079
	CD20 **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE Marked	N/A	LP-HIST-0080, LP-HIST-0079
	CD3 **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE Marked	N/A	LP-HIST-0080, LP-HIST-0079
	CD31 **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB	CE Marked	N/A	LP-HIST-0080, LP-HIST-0079

		Microscope			
CD45 **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE Marked	N/A	LP-HIST-0080, LP-HIST-0079
CD68 **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE Marked	N/A	LP-HIST-0080, LP-HIST-0079
CEA **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE Marked	N/A	LP-HIST-0080, LP-HIST-0079
Ck20 **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE Marked	N/A	LP-HIST-0080, LP-HIST-0079
Ck7 **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE Marked	N/A	LP-HIST-0080, LP-HIST-0079
CMV **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	Non-CE	N/A	LP-HIST-0080, LP-HIST-0079
Desmin **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red,	CE Marked	N/A	LP-HIST-0080, LP-HIST-0079

			optiView DAB Microscope			
	Dual Stain P16/ki67 **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE marked	N/A	LP-HIST-0080, LP- HIST-0079
	EMA **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE Marked	N/A	LP-HIST-0080, LP- HIST-0079
	GFAP **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE Marked	N/A	LP-HIST-0080, LP- HIST-0079
	Inhibin alpha **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE Marked instrument	N/A	LP-HIST-0080, LP- HIST-0079
	Ki-67 **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE Marked	N/A	LP-HIST-0080, LP- HIST-0079
	p16 **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE Marked	N/A	LP-HIST-0080, LP- HIST-0079
	p53 **2,4		Ventana Benchmark Ultra, ultraView DAB,	CE Marked	N/A	LP-HIST-0080, LP- HIST-0079

			ultraView Red, optiView DAB Microscope			
	P63 *2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE marked	N/A	LP-HIST-0080, LP- HIST-0079
	PanCk (AE1/AE3) **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE Marked	N/A	LP-HIST-0080, LP- HIST-0079
	S100 **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE Marked	N/A	LP-HIST-0080, LP- HIST-0079
	SMA *2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE Marked Instrument	N/A	LP-HIST-0080, LP- HIST-0079
	Vimentin **2,4		Ventana Benchmark Ultra, ultraView DAB, ultraView Red, optiView DAB Microscope	CE Marked	N/A	LP-HIST-0080, LP- HIST-0079
1052 Cytopathology - .02 Non gynaecological cytology	H&E and PAP **2,4	Non-gynaecological fluids Human Tissue / Cell Blocks	Manual, Thin prep 2000, Sakura Tissue-Tek VIP 6 (E2) & VIP 6 (AI) DAKEWE DP360 stainer + DAKEWE	CE Marked instrument	N/A	LP-HIST-0046 LP-HIST-0079 LP-HIST-0081 LM-HIST-0008

			CS500 coverslipper			
--	--	--	--------------------	--	--	--

Category: A

Medical pathology field - Test	Test/assay	Specimen Type	Equipment/Technique	Method (CE/Non-CE/In house developed/based on standard method)	Range of measurement	Std. ref & SOP
1012 Preparation of films on glass slides followed by microscopic examination with or without fixation and staining with dyes as required - .01 Microscopic examination for general bacteriology purposes (including enumeration and description of human cells)	Microscopic examination of General Bacteriology **1,2,3,4	All swabs, bronchial lavage / brushings, endotracheal tubes, other body fluids, pus, pleural fluids, sputum, tissue, urine, blood, CSF, middle ear, effusion, cannula, associated exudates, NICU surveillance, adverse transfusion environmental specimens, semen for culture.	Light microscope	Microscopic examination with or without fixation and staining with dyes for enumeration and description of human cells.	N/A	LM-MICRO-0001, 0003, 0004, 0005, 0006, 0007, 0008, 0010, 0011, 0012, 0014, 0016, 0019, 0022, 0023, 0033, 0034, 0035, 0050, 0054, 0055, 0061, 0063 & 0064.
1012 Preparation of films on glass slides followed by microscopic examination with or without fixation and staining with dyes as required - .03 Microscopic examination for fungi	Microscopic examination of Fungi **1,2,3,4		Light microscope	Microscopic examination with or without fixation and staining with dyes for enumeration and description of human cells.	N/A	LM-MICRO-0001, 0003, 0004, 0005, 0006, 0007, 0008, 0010, 0011, 0012, 0014, 0016, 0019, 0022, 0023, 0033, 0034, 0035, 0050, 0054, 0055, 0061, 0063 & 0064.

1013 Culture of organisms in liquid or agar based culture media with visual or instrument monitoring for growth - .01 Culture of general bacteria	Adult & Paediatric Blood Cultures*2	Blood Culture	BacTAlert	CE Marked	Growth or No Growth	LM-MICRO-0023
	Culture of General Bacteria **1,2,3,4	All swabs, bronchial lavage / brushings, endotracheal tubes, other body fluids, pus, pleural fluids, sputum, tissue, urine, blood, CSF, middle ear, effusion, cannula, associated exudates, NICU surveillance, adverse transfusion environmental specimens, semen for culture.	BacTAlert 3D System / Microscope	BacTAlert / Manual	N/A	LM-MICRO-0001, 0003, 0004, 0005, 0006, 0007, 0008, 0010, 0011, 0012, 0014, 0016, 0019, 0022, 0023, 0033, 0034, 0035, 0050, 0054, 0055, 0061, 0063 & 0064.
	Samples for Culture and Susceptibility*4	Swabs, tissue, urine, tips, sputa, blood culture, CSF, Fluids from Sterile sites, breast milk, semen, environmental samples	Culture of Organisms	CE Marked	Growth or No Growth	LM-MICRO-0015
1013 Culture of organisms in liquid or agar based culture media with visual or instrument monitoring for growth - .02 Culture of fungi	Adult & Paediatric Blood Cultures*2	Blood Culture	BacTAlert	CE Marked	Growth or No Growth	LM-MICRO-0023
	Culture of Fungi **1,2,3,4	All swabs, bronchial lavage / brushings, endotracheal tubes, other body fluids, pus, pleural fluids, sputum, tissue, urine, blood, CSF, middle ear, effusion, cannula, associated exudates, NICU surveillance, adverse transfusion environmental specimens, semen for culture.	Culture of Organisms	CE Marked	N/A	LM-MICRO-0001, 0003, 0004, 0005, 0006, 0007, 0008, 0010, 0011, 0012, 0014, 0016, 0019, 0022, 0023, 0033, 0034, 0035, 0050, 0054, 0055, 0061, 0063 & 0064.

	Samples for Culture and Susceptibility*4	Swabs, tissue, urine, tips, sputa, blood culture, CSF, Fluids from Sterile sites, breast milk, semen, environmental samples	Culture of Organisms	CE Marked	Growth or No Growth	LM-MICRO-0015
1014 Detection of bacterial, parasite, viral or fungal antigens using specific antibodies and appropriate techniques - .01 Slide agglutination,	Identification of Organisms **1,2,4	All swabs ,bronchiallavage/brushings,endotracheal tubes,other body fluids,pus,pleural fluids,sputum,tissue,urine,blood,CSF,middle ear effusion,cannual,associated exudates,NICU surveillance,adverse transfuison enviromental specimens,semen for culture	Strep grouping kit,staphaurexplus, Oxidase test, Catalase test, Nitrocefin kits	CE Marked Commercial Kits	N/A	LM-MIRO-0026,LM-MICRO-0062,LM-MICRO-0029,LM-MICRO0030,LM-MICRO-0057,LM-MICRO-0024.
1014 Detection of bacterial, parasite, viral or fungal antigens using specific antibodies and appropriate techniques - .02 Particle agglutination			Strep grouping kit,staphaurexplus, Oxidase test, Catalase test, Nitrocefin kits	CE Marked Commercial Kits	N/A	LM-MIRO-0026,LM-MICRO-0062,LM-MICRO-0029,LM-MICRO0030,LM-MICRO-0057,LM-MICRO-0024.
1014 Detection of bacterial, parasite, viral or fungal antigens using specific antibodies and appropriate techniques - .03 Enzyme immunoassay,	C. diff Quik Chek Complete	Stool	Manual	CE Marked Immunoassy Kit	Detected or Not Detected	LM-MICRO-0089
1015 Detection and/or identification of bacterial, parasite, fungal and viral nucleic acids - .01 Nucleic acid probe hybridization, CE marked commercial systems	C. diff toxin, **1,2,4	Stools	GeneXpert	PCR	N/A	LM-MICRO-00713

1015 Detection and/or identification of bacterial, parasite, fungal and viral nucleic acids - .03 Nucleic acid amplification tests, CE marked commercial systems	BV plus PCR Assay**1,2,3,4	Genital Swabs	Hamilton Seegen/ Biorad CFX96	CE marked	Normal/Positive/ Intermediate	LM-MICRO-00022,0023
	4-Plex (Flu A, Flu B, RSV and SARS-CoV-2)	Deep Nasal Swab	GeneXpert	CE	Detected or Not Detected	LM-MICRO-0072 and LM-MICRO-0086
	Biofire FilmArray BCID2 **2,3,4	Positive Blood Cultures	Biofire FilmArray	Molecular	Detected or Not Detected	LM-MICRO-0081
	Detection of CT/NG PCR using GeneXpert **1,2,3,4	Urine, genital swabs and eye swabs	GeneXpert	PCR CE Marked	Detected or Not Detected	LM-MICRO-0080
	Detection of CT/NG/TV/MG PCR using the Seegene **1,2,3,4	First void urines and Genital Swabs	Seegene Hamilton Nimbus and Biorad CFX97	PCR CE Marked	Qualitative and Quantitative	LM-SERO-0023
	Detection of GBS on the GeneXpert **1,2,3,4	LVS and Rectal Swabs	GeneXpert	PCR CE Marked	Detected or Not Detected	LM-MICRO-0075
	Molecular Detection of HSV1/2	Skin, Genital or Eye	Luminex	CE	Detected or Not Detected	LM-MICRO-0091
	Processing Positive Blood Cultures on the BioFire FilmArray **1,2,3,4	Blood Culture	BioFire FilmArray	PCR CE Marked	Detected or Not Detected	LM-MICRO-0081
	Testing for Norovirus GI and GII by PCR	Faeces	Cepheid genexpert	CE Marked	Detected / Not Detected	LM-MICRO-0084

	Testing for Trichomonas vaginalis by PCR**1,2,3,4	Genital swabs, eye swabs and urine	Cepheid genexpert	CE KIT	Detected / Not Detected	LM-MICRO-0082
	Testing of CPE isolates **1,2,3,4	Bacterial isolates	Cepheid genexpert	Ce Marked	Detected / Not Detected	LM-MICRO-0083
	Testing of respiratory samples using RP2 plus panel**1,2,3,4	NPA	Bio Fire Film Array	CE Marked	Detected / Not Detected	LM-MICRO-0083
	Xpert Xpress GBS Testing on the GeneXpert **2,3,4	Vaginal / Rectal Swabs	GeneXpert	Molecular	Detected or Not Detected	LM-MICRO-0075
1016 Identification of cultured bacteria and fungi using non-nucleic acid based techniques - .04 Identification using MALDI-TOF Spectroscopy	Identification of Organisms **1,2,4	All swabs, bronchial lavage / brushings, endotracheal tubes, other body fluids, pus, pleural fluids, sputum, tissue, urine, blood, CSF, middle ear, effusion, cannula, associated exudates, NICU surveillance, adverse transfusion environmental specimens, semen for culture.	Vitek MaldiToF MS	Spectroscopy	N/A	LM-MICRO-0076 & 0033
			Vitek 2 Compact	Biochemical testing		LM-MICRO-0058,0059,0060 & 0065. LP-MICRO-0027
1017 Measurement of antimicrobial activity and application of clinical interpretive criteria to general bacteria (rapidly growing aerobes) - .01 Anaerobes	Aerobes / Anaerobes **1,2,4		Vitek 2 Compact	Biochemical Testing	N/A	LM-MICRO-0018, 0058 & 0059
			Vitek 2 Compact and Manual	MIC / Disk Diffusion	N/A	LM-MICRO-0018, 0058 & 0059
1018 Detection of antibody response to infection using appropriate CE marked	Vz IgG	Serum	Vidas Automated	CE Marked	Immune or Non Immune	LM-SERO-0003

commercial techniques - .02 Enzyme immunoassay, using CE marked commercial systems						
1018 Detection of antibody response to infection using appropriate CE marked commercial techniques - .04 Line immunoassay, using CE marked commercial systems	Hepatitis B Core Antibodies Anti-HBc **1,2,3,4		Abbott Architect	Immunoassay	N/A	LM-SERO-0012
	Hepatitis B Surface Antigen (HBsAg) **1,2,3,4		Mini Vidas	Immunoassay	N/A	LM-SERO-0006
	Hepatitis C Antibodies Anti-HCV **1,2,3,4		Abbott Architect	Immunoassay	N/A	LM-SERO-0009
	HIV Combo Ab/Ag **1,2,3,4		Abbott Architect	Immunoassay	N/A	LM-SERO-0010
	Rubella IgG **1,2,3,4		Abbott Architect	Immunoassay	N/A	LM-SERO-0025
	Syphilis TP (Treponema pallidum) **1,2,3,4		Abbott Architect	Immunoassay	N/A	LM-SERO-0014
			Mini Vidas	Immunoassay	N/A	LM-SERO-0019
			Abbott Architect	Immunoassay	N/A	LM-SERO-0013
1029 Miscellaneous - .01 Semen	Semen Analysis **1,4	Semen	Sperminator Equipment & Manual	Count, volume, motility. CE Marked	N/A	LM-MICRO-0074

analysis (Screening Test)						
1029 Miscellaneous - .99 Miscellaneous tests	Pregnancy Testing **4	Urine	Alere Diagnostic Stick Chromatography immunoassay	CE Marked	N/A	LM-MICRO-0025

The hospital pathology laboratory is accredited for the provision of Point of Care testing in accordance with ISO 15189:2012 and ISO 22870:2017, for the tests identified in category B.

The laboratory has been awarded flexible scope in the scope classifications as noted in the scope document and in accordance with the laboratories approved and documented procedures.

Note 1 - Range may be extended for the test

Note 2 – New parameters / tests may be added

Note 3 – New matrices may be added

Note 4 – Changes to equipment / kits where the underlying methodology does not change

For further details please refer to the laboratories 'Master list of Flexible scope changes', available directly from the laboratory.